

# OptiPrep™ Application Sheet S51

## Isolation of optical system membranes and structures

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ **To access other Application Sheets** referred to in the text: return to the **2020SMemapp** file and select the appropriate S-number.

**This Application Sheet covers the purification of the following particles:**

**Section 1: Retinal rod outer segments (ROS) and ROS disks**

**Section 2: Photoreceptor outer segments**

**Section 3a: References**

**Section 3b: Recent references (after 2018)**

### 1. Retinal rod outer segments (ROS) and ROS disks

- ◆ Important technical notes, information regarding alternative methodologies and membrane analysis are contained in the “Technical Notes and Review” section (Section 1.5).
- ◆ Section 1.5 includes information on some variations to the given format and some of the more recent references

#### 1.1 Background

The use of OptiPrep™ for the isolation of mouse retinal rod outer segments (ROS) was first published in 1998 [1]. After removal of the retinas, they are repeatedly vortexed in a Ringer’s solution containing iodixanol (usually, but not always, 8% v/v/ OptiPrep™) and centrifuged at low speed. During this procedure the ROS are released into the supernatants and subsequently banded in a two-step 10%, 18% (v/v) OptiPrep™ gradient. The method was extended without modification to the isolation of ROSs from human tissue [2]. Howes et al [3] and Calvert et al [4] replaced the two-layer iodixanol gradient with a three layer one for the purification of the ROS, although the densities of the three layers were slightly different and Ringer’s solution was replaced by Locke’s solution.

More recently Liang et al [5] introduced a continuous iodixanol gradient for the purification of the ROSs and subsequently, after osmotic lysis of the ROSs, the optical disks were recovered on a second continuous iodixanol gradient. For convenience, the methodology described in this Application Sheet is adapted from ref 5 but the extraction and gradient variations for ROS purification are presented in Section 5.3.

#### 1.2 Solutions required (see Section 5.1)

- OptiPrep™
- Ringers Solution: 130 mM NaCl, 3.6 mM KCl, 2.4 mM MgCl<sub>2</sub>, 1.2 mM CaCl<sub>2</sub>, 0.02 mM EDTA, 10 mM Hepes-NaOH, pH 7.4
- Lysis buffer: 2 mM Tris-HCl, pH 7.4

#### 1.3 Centrifuge requirements

Swinging-bucket rotor for a high-speed centrifuge with a tube volume of approx 15 ml

Keep the following stock solutions at 4°C:  
 500 mM Hepes (free acid): 11.9 g per 100 ml water  
 10 mM EDTA (Na<sub>2</sub>•2H<sub>2</sub>O): 0.372 g per 100 ml water  
 1 M NaCl: 5.84 g per 100 ml of water  
 100 mM KCl: 0.745 g per 100 ml of water  
 100 mM MgCl<sub>2</sub>•6H<sub>2</sub>O: 2.03 g per 100 ml of water  
 100 mM CaCl<sub>2</sub>•2H<sub>2</sub>O: 1.47 g per 100 ml water  
 100 mM Tris (free base): 1.21 g per 100 ml

Solution B: Mix together the following volumes of the above stocks, **except the Tris** (in order), 4 ml, 0.4 ml, 26 ml, 7.2 ml, 4.8 ml and 2.4 ml; adjust to pH 7.4 with NaOH and make up to 200 ml

Solution C: Add 4 ml of Tris stock solution to 100 ml of water; adjust to pH 7.4 with 1 M HCl; make up to 200 ml

#### 1.4 Protocol

1. Make an 8% (v/v) OptiPrep™ solution by diluting 8 ml of Solution A to 100 ml with Solution B.

2. Suspend 12 mouse retinas in 0.12 ml of this solution. Carry out this and the following steps at room temperature.
3. Vortex for 1 min and then centrifuge at 200 g for 1 min.
4. Carefully recover the supernatant (without disturbing the pellet) and retain.
5. Resuspend the pellet in 0.12 ml of the 8% OptiPrep™ solution and repeat steps 3-4.
6. Repeat this procedure five or six times.
7. Prepare two gradient solutions of 10% and 30% (v/v) OptiPrep™ by diluting 1 ml and 3 ml of Solution A to 10 ml with Solution B respectively.
8. Using a two-chamber gradient maker or a Gradient Master™ prepare a 10-12 ml continuous gradient from these two solutions in tubes for the swinging-bucket rotor ([see Section 5.2](#)).
9. Layer the combined supernatants (approx. 1.5 ml) on top and centrifuge at 26,500 g for 30 min.
10. Collect the band of ROSs (about two thirds of the way from the top of the gradient).
11. Dilute with 3 vol. of Solution B and centrifuge at 500 g for 3 min to pellet the nuclei.
12. Harvest the ROSs by centrifuging at 26,500 g for 30 min.
13. To prepare the optical disks, lyse the ROSs by suspending them in 2 ml of Solution C and maintain them at 0°C for 15 h.
14. Prepare two gradient solutions of 15% and 40% (v/v) OptiPrep™ from 1.5 ml and 4 ml of Solution A, each diluted to 10 ml with Solution B respectively.
15. Using a two-chamber gradient maker or a Gradient Master™ prepare a 10-12 ml continuous gradient from these two solutions in tubes for the swinging-bucket rotor ([see Section 5.2](#)).
16. Layer the lysed ROSs on top of the gradient and centrifuge at 26,500 g for 30 min.
17. Harvest the disks that band about two thirds of the way from the top of the gradient

## 1.5 Technical Notes and Review

### 1.5.1 Homogenization media and gradient solutions

All iodixanol solutions in this Optiprep™ Application Sheet are prepared simply by diluting Optiprep™ with Ringers Solution, which will mean that the levels of KCl, MgCl<sub>2</sub>, CaCl<sub>2</sub>, EDTA and buffer will be reduced. If it is considered important that these levels should be maintained in all density solutions then a working solution of 50% (w/v) iodixanol should be first prepared by diluting 5 vol. of Optiprep™ with 1 vol. of 130 mM NaCl, 21.6 mM KCl, 14.4 mM MgCl<sub>2</sub>, 7.2 mM CaCl<sub>2</sub>, 0.12 mM EDTA, 60 mM Hepes-NaOH, pH 7.4. The working solution would thus contain the same concentrations of KCl, MgCl<sub>2</sub>, CaCl<sub>2</sub>, EDTA and Hepes as the Ringers buffer. Further dilutions with Ringers buffer should then need to be adjusted appropriately from those given in the protocol to give gradient solutions of the correct density. Note that if the concentration of NaCl in the diluent were also raised then the osmolality of the 50% (w/v) iodixanol will be unacceptably high. Strategies for preparing working solutions for mammalian tissues and cells are given in [Application Sheet S01](#).

### 1.5.2 Construction of gradients

If neither of these gradient-making devices is available then a continuous gradient can be prepared by diffusion of a discontinuous gradient. For example, the 10-30% (v/v) Optiprep™ gradient (step 7) might be generated from a discontinuous gradient of equal volumes of 10%, 18%, 24% and 30% Optiprep™. Likewise the 15-40% (v/v) Optiprep™ gradient (step 14) might be generated from a discontinuous gradient of equal volumes of 15%, 21%, 27%, 33% and 40% Optiprep™. For more information on gradient construction see [Application Sheet S03](#).

### 1.5.3 Variations in gradient and centrifugation format

#### Continuous gradients

Continuous iodixanol gradients normally conform to the above protocol [5-13], although Peshenko et al [14,15] used a 20-30% (v/v) OptiPrep™ gradient and a higher g-force of 75,000 g, for 40 min.

#### Discontinuous gradients

The commonly used discontinuous gradient format involves layering the extracted ROSs in 8% (v/v) Optiprep™ over a two layers of 10% and 18% (v/v) Optiprep™; the centrifugation conditions are however somewhat diverse: 3,300 g for 10 min [1,2], 50,000 g for 1 h [16] 70,000 g for 1 h [17]. Variants in the density of the gradient layers include 8.3% and 20% with the sample in 6.7%, centrifuged at only 1,425 g, for 15 min [4]. Gradients comprising three layers: 10, 20 and 30% (sample in 8%), after centrifugation at 17,000 g for 50 min the ROSs banded at the 20%/30% interface [3] and four layers: 8%, 12%, 16% and 20% (sample in 2%) [18,19].

- ◆ Occasionally after repeated extractions with 8% (v/v) Optiprep™ no gradient is used [20].
- ◆ Refs 21-23 are meetings abstracts, reporting the use of Optiprep™ but no details are provided.
- ◆ **Refs 24-35 all report the use of iodixanol gradients prepared from OptiPrep™ for ROS isolation using methods identical or very similar to those described above**

## 2. Photoreceptor outer segments (POS)

Jiang et al [36] and Hazim et al [37] described homogenizing retinas in 130 mM NaCl, 3.6 mM KCl, 2.4 mM MgCl<sub>2</sub>, 1.2 mM CaCl<sub>2</sub>, 0.02 mM EDTA, and 10 mM HEPES, pH 7.4; then after centrifugation at 100 g for 1 min, the supernatant was layered over 8, 10 and 15% (w/v) iodixanol and centrifuged at 12,000 g for 20 min. The outer segments banded at the 10-15% interface. In an interesting variation Pelkonen et al [38] described homogenizing the tissue in 8% iodixanol (in a bicarbonate buffer). After an initial centrifugation at 720 g for 3 min (repeated five times); the suspension was diluted to 2% iodixanol and layered over 10% and 20% iodixanol, banding the POS at 14,000 g for 30 min. Rao et al [39] used a similar approach to that described in ref 38 but the gradient was a continuous one (10-30% iodixanol).

### 3a. References

1. Tsang, S.H., Burns, M.E., Calvert, P.D., Gouras, P., Baylor, D.A., Goff, S.P. and Arshavsky, V.Y. (1998) *Role for the target enzyme in deactivation of photoreceptor G protein in vivo* Science, **282**, 117-121
2. Mata, N.L., Weng, J. and Travis, G.H. (2000) *Biosynthesis of a major lipofuscin fluorophore in mice and humans with ABCR-mediated retinal and macular degeneration* Proc. Natl. Acad. Sci. USA, **97**, 7154-7159
3. Howes, K.A., Pennesi, M.E., Sokal, I., Church-Kopish, J., Schmidt, B., Margolis, D., Frederick, J.M., Rieke, F., Palczewski, K., Wu, S.M., Detwiler, P.B. and Baehr, P. (2002) *GCAP1 and rescues rod photoreceptor response in GCAP1/GCAP2 knockout mice* The EMBO J., **21**, 1545-1554
4. Calvert, P.D., Govardovskii, V.I., Krasnoperova, N., Anderson, R.E., Lem, J. and Makino, C.L. (2001) *Membrane protein diffusion sets the speed of rod phototransduction* Nature, **411**, 90-94
5. Liang, Y., Fotiadis, D., Filipek, S., Saperstein, D. A., Palczewski, K. and Engel, A. (2003) *Organization of the G protein-coupled receptors rhodopsin and opsin in native membranes* J. Biol. Chem., **278**, 21655-21662
6. Fotiadis, D., Liang, Y., Filipek, S., Saperstein, D. A., Engel, A. and Palczewski, K. (2004) *The G protein-coupled receptor rhodopsin in the native membrane* FEBS Lett., **564**, 281-288
7. Liang, Y., Fotiadis, D., Filipek, S., Saperstein, D.A., Engel, A. and Palczewski, K. (2004) *Rhodopsin signaling and organization in heterozygote rhodopsin knockout mice* J. Biol. Chem., **279**, 48189-48196

8. Maeda, A., Maeda, T., Imanishi, Y., Kuksa, V., Alekseev, A., Bronson, J.D., Zhang, H., Zhu, L., Sun, W., Saperstein, D.A., Ricke, F., Baehr, W. and Palczewski, K. (2005) *Role of photoreceptor-specific retinal dehydrogenase in the retinoid cycle in vivo* J. Biol. Chem., **280**, 18822-18832
9. Saperstein, D.A., Fotiadis, D., Liang, Y., Filipek, S., Palczewski, K. and Engel, A. (2003) *The structure of murine outer segment disk membranes using atomic force microscopy* Invest. Ophthalmol Vis. Sci., **44**, E-abstract, 3175
10. Nickell, S., Park, P.S-H., Baumeister, W. and Palczewski, K. (2007) *Three-dimensional architecture of murine outer rod segments determined by cryoelectron tomography* J. Cell Biol., **177**, 917-925
11. Mustafi, D., Kevany, B.M., Genoud, C., Okano, K., Cideciyan, A.V., Sumaroka, A., Roman, A.J., Jacobson, S.G., Engel, A., Adams, M.D. and Palczewski, K. (2011) *Defective photoreceptor phagocytosis in a mouse model of enhanced S-cone syndrome causes progressive retinal degeneration* FASEB J. **25**, 3157–3176
12. Ziccardi, L., Vijayarathy, C., Bush, R.A. and Sieving, P.A. (2012) *Loss of retinoschisin (RS1) cell surface protein in maturing mouse rod photoreceptors elevates the luminance threshold for light-driven translocation of transducin but not arrestin* J. Neurosci., **32**, 13010–13021
13. Gilliam, J.C., Chang, J.T., Sandoval, I.M., Zhang, Y., Li, T., Pittler, S.J., Chiu, W. and Wensel, T.G. (2012) *Three-dimensional architecture of the rod sensory cilium and its disruption in retinal neurodegeneration* Cell. **151**, 1029–1041
14. Peshenko, I.V., Olshevskaya, E.V., Savchenko, A.B., Karan, S., Palczewski, K., Baehr, W. and Dizhoor, A.M. (2011) *Enzymatic properties and regulation of the native isozymes of retinal membrane guanylyl cyclase (RetGC) from mouse photoreceptors* Biochemistry, **50**, 5590–5600
15. Peshenko, I.V., Olshevskaya, E.V., Azadi, S., Molday, L.L., Molday, R.S. and Dizhoor, A.M. (2011) *Retinal degeneration 3 (RD3) protein inhibits catalytic activity of retinal membrane guanylyl cyclase (RetGC) and its stimulation by activating proteins* Biochemistry, **50**, 9511– 9519
16. Dizhoor, A.M. Woodruff, M.L., Olshevskaya, E.V., Cilluffo, M.C., Cornwall, M.C., Sieving, P.A. and Fain, G.L. (2008) *Night blindness and the mechanism of constitutive signaling of mutant G90D rhodopsin* J. Neurosci., **28**, 11662-11672
17. Burns, M.E., Mendez, A., Chen, C-K., Almuete, A., Quillinan, N., Simon, M.I., Baylor, D.A. and Chen, J. (2006) *Deactivation of phosphorylated and nonphosphorylated rhodopsin by arrestin splice variants* J. Neurosci., **26**, 1036-1044
18. Krispel, C.M., Chen, D., Melling, N., Chen, Y-J., Martemyanov, K.A., Quillinan, N., Arshavsky, V.Y., Wensel, T.G., Chen, C-K. and Burns, M.E. (2006) *RGS expression rate-limits recovery of rod photoresponses* Neuron, **51**, 409-416
19. Gross, O.P. and Burns, M.E. (2010) *Control of rhodopsin's active lifetime by arrestin-1 expression in mammalian rods* J. Neurosci., **30**, 3450–3457
20. Nair, K.S., Hanson, S.M., Mendez, A., Gurevich, E.V., Kennedy, M.J., Shestopalov, V.I., Vishnivetskiy, S.A., Chen, J., Hurley, J.B., Gurevich, V.V. and Slepak, V.Z. (2005) *Light-dependent redistribution of arrestin in vertebrate rods is an energy-independent process governed by protein-protein interactions* Neuron, **46**, 555-567
21. Wang, Q., Hu, G., Leitges, M. and Wensel, T.G. (2006) *Phosphorylation of phototransduction GAP RGS9-1 depends on the  $\alpha$  isoform of protein kinase C* Invest. Ophthalmol. Vis. Sci., **47**, E-Abstr 822
22. Finnemann, S.C. and Chang, Y. (2006) *Regulation of RPE phagocytosis by integrin receptor-tetraspanin surface membrane domains* Invest. Ophthalmol. Vis. Sci., **47**, E-Abstr 822
23. Coleman, J.A., Djajadi, H.R., Molday, L.L. and Molday, R.S. (2012) *Disruption of the P<sub>r</sub>-ATPase aminophospholipid flippase Atp8a2 gene suggests a role for phosphatidylserine in photoreceptor outer segments* Invest. Ophthalmol. Vis. Sci., **53**, Abstr. 743- 620
24. Moaven, H., Koike, Y., Jao, C.C., Gurevich, V.V., Langen, R. and Chen, J. (2013) *Visual arrestin interaction with clathrin adaptor AP-2 regulates photoreceptor survival in the vertebrate retina* Proc. Natl. Acad. Sci. USA, **110**, 9463–9468
25. Skiba, N.P., Spencer, W.J., Salinas, R.Y., Lieu, E.C., Thompson, J.W. and Arshavsky, V.Y. (2013) *Proteomic identification of unique photoreceptor disc components reveals the presence of PRCD, a protein linked to retinal degeneration* J. Proteome Res., **12**, 3010-3018
26. Coleman, J.A., Zhu, X., Djajadi, H.R., Molday, L.L., Smith, R.S., Libby, R.T., John, S.W.M. and Molday, R.S. (2014) *Phospholipid flippase ATP8A2 is required for normal visual and auditory function and photoreceptor and spiral ganglion cell survival* J. Cell Sci., **127**, 1138–1149
27. Wensel, T.G. and Gilliam, J.C. (2015) *Three-dimensional architecture of murine rod cilium revealed by cryo-EM* In Methods Mol. Biol., **1271**, Rhodopsin: Methods and Protocols (ed. Jastrzebska, B.), Springer Science+Business Media, New York, pp 267-292
28. Rakshit, T., Senapati, S., Sinha, S., Whited, A.M. and Park, P.S-H. (2015) *Rhodopsin forms nanodomains in rod outer segment disc membranes of the cold-blooded Xenopus laevis* PLoS One, **10**: e0141114

29. Calzia, D., Panfoli, I., Heinig, N., Schumann, U., Ader, M., Traverso, C.E., Funk, R.H.W. and Roehlecke, C. (2016) *Impairment of extramitochondrial oxidative phosphorylation in mouse rod outer segments by blue light irradiation* Biochimie, **125**, 171-178
30. He, F., Agosto, M.A., Anastassov, I.A., Tse, D.Y., Wu, S.M. and Wensel, T.G. (2016) *Phosphatidylinositol-3-phosphate is light-regulated and essential for survival in retinal rods* Sci. Rep., **6**: 26978
31. Berry, J., Frederiksen, R., Yao, Y., Nymark, S., Chen, J. and Cornwall, C. (2016) *Effect of rhodopsin phosphorylation on dark adaptation in mouse rods* J. Neurosci., **36**, 6973–6987
32. Peshenko, I.V., Olshevskaya, E.V. and Dizhoor, A.M. (2016) *Functional study and mapping sites for interaction with the target enzyme in retinal degeneration 3 (RD3) protein* J. Biol. Chem., **291**, 19713–19723
33. Rajala, R.V.S., Rajala, A., Kooker, C., Wang, Y. and Anderson, R.E. (2016) *The Warburg effect mediator pyruvate kinase M2 expression and regulation in the retina* Sci. Rep., **6**: 37727
34. Rakshit, T., Senapati, S., Parmar, V.M., Sahu, B., Maeda, A. and Park, P.S-H. (2017) *Adaptations in rod outer segment disc membranes in response to environmental lighting conditions* BBA – Mol. Cell Res., **1864**, 1691-1702
35. Maity, S., Ilieva, N., Laio, A., Torre, V. and Mazzolini, M. (2017) *New views on phototransduction from atomic force microscopy and single molecule force spectroscopy on native rods* Sci. Rep., **7**: 12000
36. Jiang, M., Esteve-Rudd, J., Lopes, V.S., Diemer, T., Lillo, C., Rump, A. and Williams, D.S. (2015) *Microtubule motors transport phagosomes in the RPE, and lack of KLC1 leads to AMD-like pathogenesis* J. Cell Biol., **210**, 595–611
37. Hazim, R., Jiang, M., Esteve-Rudd, J., Diemer, T., Lopes, V.S. and Williams, D.S. (2015) *Live-cell imaging of phagosome motility in primary mouse RPE cells* In Retinal Degenerative Diseases, Adv. Exp. Med. Biol., **854**, (ed. C. Bowes Rickman et al.) Springer International Publishing Switzerland, pp 751-755
38. Pelkonen, L., Sato, K., Reinisalo, M., Kidron, H., Tachikawa, M., Watanabe, M., Uchida, Y., Urtti, A and Terasaki, T. (2017) *LC-MS/MS based quantitation of ABC and SLC transporter proteins in plasma membranes of cultured primary human retinal pigment epithelium cells and immortalized ARPE19 cell line* Mol. Pharmaceutics, **14**, 605-613

### 3b. Recent references

1. Makino, C.L., Duda, T., Pertzev, A. and Sharma, R.K. (2018) *Experimental approaches for defining the role of the Ca<sup>2+</sup>-modulated ROS-GC system in retinal rods of mouse* In Mouse Retinal Phenotyping: Methods and Protocols, Methods in Mol. Biol., **1753** (ed. Tanimoto, N.), Springer Science+Business Media, LLC, pp 129-158
2. Senapati, S. and Park, P.S-H. (2019) *Investigating the nanodomain organization of rhodopsin in native membranes by atomic force microscopy* In Atomic Force Microscopy: Methods and Protocols, Methods in Molecular Biology, vol. **1886** (ed Santos, N.C. and Carvalho, F.A. Springer Science+Business Media LLC New York, pp 61-74
3. Dilan, T.L., Moye, A.R., Salido, E.M., Saravanan, T., Kolandaivelu, S., Goldberg, A.F.X. and Ramamurthy, V. (2019) *ARL13B, a Joubert syndrome-associated protein, is critical for retinogenesis and elaboration of mouse photoreceptor outer segments* J. Neurosci., **39**, 1347–1364
4. Spencer, W.J., Ding, J-D., Lewis, T.R., Yu, C., Phan, S., Pearing, J.N., Kim, K-Y., Thor, A., Mathew, R. et al (2019) *PRCD is essential for high-fidelity photoreceptor disc formation* Proc. Natl. Acad. Sci., **116**, 13087–13096

**OptiPrep™ Application Sheet S51; 8<sup>th</sup> edition, January 2020**