

OptiPrep™ Application Sheet M09

Analysis of protein complex formation, microtubules and cytoskeleton

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ To access other Application Sheets referred to in the text: return to the 2020Macroapp file and select the appropriate M-number.

1. Background

Determination of the sedimentation coefficient of a protein, traditionally in a 5-20% sucrose gradient, in order to determine its molecular mass, is a technique that was developed many years ago. Although alternative methods, notably polyacrylamide gel electrophoresis, have achieved a wide popularity, there are certain situations that are better suited to density gradients, for example the analysis of heavily glycosylated proteins, which run anomalously on gels. Another area that may be better suited to density gradient analysis is the study of protein-protein interactions.

The standard sucrose gradient analysis, which is carried out in swinging-bucket rotors for 4-16 h, may be less than ideal for determining protein-protein interactions. Macromolecular complexes may be insufficiently stable to survive these long centrifugation times and they may be intolerant of the high hydrostatic pressures generated in a swinging-bucket rotor [1,2]. In addition Timasheff [3] pointed out that solutions of low water activity (high osmolality) can remove bound water from proteins and could cause changes in stability of the protein and its propensity to aggregate. The use Nycodenz® will reduce the osmolality (raise the water activity) of gradients considerably, while with iodixanol gradients can be made isoosmotic throughout the entire useful density range. Gradient made from one or other of these solutes may therefore be an important advantage in studying protein-protein interactions in density gradients. Some reports also highlight the functional problems associated with sucrose gradients; for example these caused proteolysis of a kinesin-related motor protein, while in iodixanol gradients there was no proteolysis whatsoever [4]. This Application Sheet covers the following topics:

- ◆ Section 2: Alzheimer's A β peptide and other neurodegenerative disease-associated proteins
- ◆ Section 3: Non-muscle myosin II (NMII)
- ◆ Section 4: Dimerization of the kinesin-related motor protein KIF1A
- ◆ Section 5: α -Synuclein aggregation
- ◆ Section 6: Microtubule fragment size analysis
- ◆ Section 7: Cytoskeleton analysis
- ◆ Section 8: Prion proteins and fibrils
- ◆ Section 9: Association of *Drosophila* proteins with lipophorin
- ◆ Section 10: Hepatitis and Herpes virus proteins
- ◆ Section 11: Other recent applications

2. Alzheimer's disease amyloid-beta (A β) peptide

2a. Aggregation (oligomerization) studies

2a-1. Introduction

Moir et al [5] studied the aggregation of the A β peptide by Zn²⁺, Cu²⁺, EDTA or a pH 5.5 buffer. The incubations were simply layered over a series of Nycodenz® barriers of increasing density; the concentrations were 43%, 44%, 44.6%, 46.8% and 47% (w/v) in microcentrifuge tubes and centrifuged at 16,000 g for 10 min. Tubes were then analyzed by rapid freezing to -170°C for 1 h and "fractionated" by cutting the tube. All the treatments, except EDTA, allowed the aggregates to sediment through the 43 and 44% barriers but only the Zn²⁺ induced aggregates sedimented through the 44.6 and 46.8% barriers.

The analysis of oligomerization of the β -amyloid A β peptide in iodixanol gradients, using a high performance near-vertical rotor, developed by Ward et al [6], may be applicable to any protein oligomerization studies and is thus given here in some detail. The short sedimentation path length of the rotor (18 mm) means that the centrifugation time can be reduced to 3 h and the hydrostatic pressure is correspondingly low. The authors used the gradients to analyze the interactions of the β -amyloid (A β) peptide that is implicated in Alzheimer's disease, but the strategy could be used or adapted to any protein oligomerization or protein-protein interaction study. The protocol described below is adapted from ref 6. **For Notes see Section 2a-5**

2a-2. Solutions required

- A. OptiPrep™
- B. OptiPrep™ diluent: for the study of the oligomerization of β -amyloid (A β) peptide [6] the diluent used was phosphate-buffered saline (PBS), but any diluent compatible with the protein-protein interaction may be used (see Note 1). **Include protease inhibitors in solutions as required.**

2a-3. Ultracentrifuge rotor requirements

A near-vertical (e.g. Beckman NVT100, NVT90 or NVT65.2) or vertical (e.g. Beckman VTi90, VTI65.2 or Sorvall TV1665 or 70V6) rotor: all these rotors accommodate tubes for approx 5 ml total volume (see Note 2).

2a-4. Protocol

1. Carry out all operations at 0-4°C or any temperature compatible with the interactions under investigation.
2. Dilute OptiPrep™ with PBS to produce the solutions used in the next step (see Note 1).
3. Layer sequentially in tubes for the chosen rotor: 0.65 ml of 50%, 40% and 30% (w/v) iodixanol, 1.95 ml of 20%, 0.65 ml of 10% and 0.3 ml of 5% (w/v) iodixanol (see Note 3)
4. Layer 0.3 ml of the protein solution on top and centrifuge at 350,000 g for 3 h at 0-4°C, using controlled acceleration and deceleration over the 0-2000 rpm range (see Note 4).
5. Collect the gradient in 0.3-0.35 ml fractions either by upward displacement with Maxidens™ and an Axis-Shield Gradient Collector, tube puncture or by automatic aspiration from the meniscus using a Labconco Auto Densi-flow® (see Notes 5 and 6).

◆ See Notes 7-11 for other examples of the use of iodixanol gradients for A β peptide studies

2a-5. Notes

1. If it is important to maintain a constant background of buffer concentration (e.g. 20 mM Tris-HCl, pH 8.0) and low concentration of some other additive (e.g. 1 mM EDTA and 1 mM DTT); then a 50% (w/v) iodixanol working solution (WS) should first be prepared by mixing 5 vol. of OptiPrep™ with 1 vol. of 6 mM EDTA, 6 mM DTT, 120 mM Tris-HCl, pH 8.0. Further dilutions are then obtained by dilution of the WS with 20 mM Tris-HCl, pH 8.0, 1 mM DTT, 1 mM EDTA. All gradient solutions will then contain the requisite concentration of buffer and additive. For more information on the preparation of density gradient solutions **see Application Sheet M01**.
2. Vertical or near-vertical rotors are ideal for sedimentation velocity separations because, (a) the sample layer achieves a very narrow zone-width after reorientation and (b) the sedimentation path length is short. If other rotors (fixed-angle or swinging-bucket) are used the sample volume must never be greater than 10% of the gradient volume and the centrifugation time will need increasing to take account of the longer sedimentation path length (particularly in larger volume swinging-bucket rotors).

3. Because of the short path length, the gradient will become more or less continuous within 30 min. An alternative approach may be to produce a linear continuous gradient in the centrifuge tube before layering the sample (see [Application Sheet M02](#)); however the irregular density profile achieved by using the non-uniform layer volumes described in Step 3 of the Protocol will be lost. The consequences of this are not clear.
4. If the method is being used to study the oligomerization of other proteins the RCF or centrifugation time may need modulation according to the size of the protein particles.
5. An example of the resolution, which is achievable with this gradient system, is given in Figure 1. Freshly-prepared hexafluoro-2-propanol (HFIP) treated β -amyloid ($A\beta$) peptide oligomerizes during incubation at 35°C and the distribution of this peptide across the gradient was analyzed by SDS-PAGE after 0 min, 30 min, 18 h and 18 days of incubation at this temperature. The position of the peptide in the gradient reflects the extent and type of association of the peptide monomers.
6. For more information on methods of harvesting gradients see [Application Sheet M04](#).
7. The method has also been adapted to small volume swinging-bucket rotors: Rzepecki et al [7] used a similar gradient, which was scaled down to use in a Beckman TLS55 (2.2 ml tubes) with centrifugation at 259,000 g for 4 h and Lockhart et al [8] used a Beckman MLS50 (5 ml tubes) for studying the formation of β -amyloid fibrils at 268,000 g for 3 h. The Rzepecki et al [7] gradient was also used by Funke et al [9] to study the effect of the D-enantiomeric peptide D3, which caused a shift of the $A\beta$ peptide from an oligomeric form to much larger aggregates and the $A\beta$ binding peptide L3 [10].
8. More recent studies of the oligomerization procedure [11-13] have also used the TLS55 rotor under more or less similar centrifugation conditions; the gradients were however slightly modified; a discontinuous gradient of 0.26 ml each of 50%, 40% and 30% iodixanol, 0.78 ml of 20%, 0.26 ml of 10% and 0.1 ml of 5% iodixanol. The gradients were used to analyse $A\beta$ (1-42) aggregates using 100 μ l of samples, which was layered on top. The gradient system described by Brener et al [11] has been used to show the preferential binding of monomeric $A\beta$ to D-peptides [13-16] and GM1 [17].
9. Sehlin et al [18] used an iodixanol gradient similar to that of Ward et al [6] each step was 0.65 ml except for the 30% (w/v) iodixanol (1.95 ml) and omitted the 5% layer. The gradient effectively resolved $A\beta$ peptide aggregates of different sizes.
10. $A\beta$ prions have also been extensively purified in a multi-step method. In the primary fractionation of a brain homogenate the latter was adjusted to 18% (w/v) iodixanol and layered over 30% and 35% iodixanol. After centrifugation at 60,000 g for 20 min the top lipid layer was discarded and the material within the denser two layers collected. After dilution with buffer, this was layered over a second gradient of 26% and 35% iodixanol and centrifuged at the same speed for 40 min. Again the top layer was discarded and the material in the two denser layers recovered prior to further processing [19].
11. More recently a simple discontinuous flotation gradient (36%, 24% and 0% iodixanol) centrifuged at 54,000 g for 3 h was used to separate less dense LDL-bound $A\beta$ from $A\beta$ [20].
12. Iodixanol gradients have recently been used to purify metal nanoparticle-conjugated $A\beta$ -specific ligands to increase binding affinity [21].




Figure 1: Iodixanol gradient sedimentation velocity analysis of β -amyloid ($A\beta$) peptide oligomerization: effect of time of incubation (low density on left). For more details see text. Reproduced from ref 6 with kind permission of the authors, ©The Biochemical Society

2b. Other neurodegenerative disease-related protein studies

2b-1. Fibrillar structures

Levy et al [22] adapted the method described in Section 2 to an analysis of islet amyloid fibril formation using a 3 ml pre-formed continuous 0-60% (w/v) iodixanol gradient, generated from equal volumes of OptiPrep™ and water (pH 2) in tubes for a Beckman SW60Ti swinging-bucket rotor. The gradients were loaded with 0.6 ml of sample and centrifuged at 150,000 g for 24 h. The gradient separated a series of low molecular mass pre-fibrillar assemblies in the top half of the gradient from the fibrils at the bottom. More recently Frenzel et al [23] separated fibrils from smaller aggregates and monomers by iodixanol density gradient centrifugation. The gradient comprised layers of 50%, 40%, 30%, 20%, 10% and 5% (v/v) OptiPrep™ iodixanol in 10 mM phosphate buffer pH 7.4; the sample was top-loaded. Centrifugation in a small volume (2.2 ml tubes) swinging-bucket rotor for 3 h at 259,000 g, resolved the three forms of the protein.

◆ See Section 3 for a related method for separating monomeric and filamentous protein forms.

The method developed by Ward et al [6] has been adapted to the slightly larger tubes (6 ml) of the Sorvall TV865 by Khlistunova et al [24] for studying the Tau protein; usually a soluble protein, but in some neurodegenerative diseases it forms “paired helical filaments” and associates with microtubules. The gradient was generated from 0.85 ml of 50% and 40%, 2.2 ml of 30%, 0.85 ml of 10% and 0.3 ml of 5% (w/v) iodixanol. Sarcosyl-insoluble fractions from neuronal tissue have also been purified in an iodixanol gradient [25]. Gradients covering approx. the same density range (145,000 g for 12 h) have analyzed the association between Tau and a 14-3-3 protein [26].

3. Non-muscle myosin II (NMII)

This investigation was reported by Shutova et al [27]. Detergent-soluble cell lysates was layered on to a discontinuous gradient of 50%, 25%, 12%, and 6% OptiPrep™ layers. After centrifugation at 80,000 rpm for 1 h the gradient resolved monomeric and filamentous forms of NMII which had sedimentation coefficients of approx 7S and 16S respectively. The monomeric form was found to be increased to a greater extent in blebbistatin-treated cells. The same gradient was reported in ref 28.

4. Kinesin-related motor protein KIF1A dimerization

Rashid et al [4] used a 3.8 ml continuous gradient of 5-40% (w/v) iodixanol overlaid by 0.2 ml of sample in a Beckman SW60Ti swinging-bucket rotor at 100,000 g_{av} for 16 h at 4°C, to identify dimers of KIF1A.

5. α -Synuclein aggregation

Rather larger protein aggregates of cytoplasmic α -synuclein from COS7 cells have been analyzed in discontinuous gradients of 2.5%, 25% and 35% (w/v) iodixanol, centrifuged at only 50,000 g for 30 min [29].

6. Microtubule fragment size analysis

MacCormick et al [30] characterized a large protein complex containing a nerve-growth-factor-activated ERK (extracellular-signal-regulated kinase) and MEK. Part of the fractionation process included the isolation of a non-ionic detergent insoluble fraction that also contained microtubule fragments that contained some bound kinases. Glycerol gradients were not effective in resolving the larger microtubule fragments; the following iodixanol gradient method devised by MacCormick et al [30] however provided excellent resolution of these particles.

1. Permeabilize PC12 cells by a single passage through a ball-bearing homogenizer.
2. Release the intracellular material by the permeabilization and separate from the residual “cells” by centrifugation at 1000 g for 10 min.
3. Concentrate the supernatant material by pelleting at 100,000 g through a cushion of 10% (w/v) sucrose in 20 mM MOPS (pH 7.2) containing 1 mM EGTA and 1 mM Na_3VO_4 .

4. Suspend the pellet in buffer containing potassium aspartate, potassium gluconate and potassium glutamate (all 38 mM), 20 mM MOPS, 10 mM potassium bicarbonate, 0.5 mM magnesium carbonate, 1 mM EDTA, 1 mM EGTA, pH 7.1.
 5. Prepare a 30% (w/v) iodixanol solution in the same buffer and make a continuous 0-30% iodixanol prepared using a two-chamber gradient maker or Gradient Master™. If neither of these devices is available make up a discontinuous gradient from equal volumes of 0%, 10%, 20% and 30% (w/v) and allow to diffuse overnight at 4°C. For more information about preparing continuous iodixanol gradients [see Application sheet M02](#).
 6. Layer the sample on top of the gradient (sample volume <10% of total gradient volume). The gradient separation is a sedimentation velocity one, so its resolving power is inversely proportional to the sample volume. Centrifuge at 200,000 g for 1-3 h (turn off the brake during deceleration from 2000 rpm or use a controlled deceleration program).
 7. Collect the gradient either by tube puncture or collection from the meniscus. For more information about gradient fractionation [see Application Sheet M05](#). MacCormick et al [30] collected the gradient in 27 fractions.
- ◆ After a centrifugation for 1 h the microtubule fragments from PC12 cells has been resolved on the basis of size into four discrete fractions, as determined by the α -tubulin distribution [30] but ERK was only associated with three of them.

7. Cytoskeleton analysis

Plasma membrane and cytoskeleton in *Drosophila* extracts have been effectively separated on 10-40% (w/v) iodixanol gradients prepared by dilution of OptiPrep™ with 0.25 M Sucrose, 10 mM Tris-HCl, pH 8.0 centrifuged at 250,000 g for 3 h [31]. It provided a very clear separation of the syntaxin 1A (plasma membrane marker) the denser myosin II heavy chain (cytoskeleton marker).

A widely used discontinuous sedimentation velocity gradient originally devised by Majoul et al [32], comprising 2.5%, 5%, 7.5%, 10%, 12.5%, 15%, 17.5%, 20%, and 30% (w/v) iodixanol was used by Chen et al [33] who studied the cytoskeleton structure in a 65,000 g (1 h) fraction from a post-nuclear supernatant of a mouse spinal cord homogenate. The gradients were allowed to diffuse at 4°C for 3-4 h solution before loading the fraction and centrifuging at 170,000 g for 90 min. Both actin and tubulin showed a biphasic distribution, being concentrated in the top quarter of the gradient but also demonstrating a sharp band about three-quarters of the way down the gradient. Mice deficient in a presenilin binding protein called “modifier of cell adhesion” (MOCA) lacked the sharp denser band.

Cofilin-actin rods have been isolated from cultured nerve cells in a double gradient strategy [34]. All steps were performed at 4°C. A low-speed supernatant was layered over 10% and 15% (w/v) iodixanol and centrifuged at 6,650 g for 10 min to band the rods above the interface between the two iodixanol layers. The collected fraction was diluted with 13% iodixanol (vol. ratio of approx 0.4:1.7) and centrifuged at 166,000 g for 2 h to create a self-generated gradient – a small volume vertical rotor is best suited to self-generated gradient formation ([see OptiPrep Application Sheet M03 for more information on self-generated gradients](#)). The rods band close to the bottom of the gradient. Once recovered, the rods in the fraction (after dilution with an equal volume of buffer) may be concentrated on a 25% iodixanol cushion at 5,000 g for 15 min [33]. Actin rods have also been analyzed on a two-layer gradient (10% and 15% v/v OptiPrep™, centrifuged at 6,600 g for 10 min [35]).

Phosphorylation of the vasodilator-stimulated phosphoprotein (VASP) controls its influence on F-actin-related processes and Lin et al [36] studied them in *Dictyostelium*. A discontinuous iodixanol gradient of 5%, 10%, 15% and 20% (w/v) iodixanol centrifuged at approx 150,000 g for 18 h. WASP banded broadly in the bottom half of the gradient but in the presence of VASP, the WASP distribution shifted significantly to a higher density.

8. Prion proteins

Tixador et al [37] solubilized the proteins from a mouse brain homogenate in a buffered solution of dodecyl- β -D-maltoside and *N*-lauryl sarcosine. The material was layered on top of a continuous 10–25% (w/v) iodixanol gradient (4.8 ml total volume) and centrifuged 285,000 g for 45 min. The host-encoded prion protein (PrP^C) banded close to the top of the gradient, while the PrP^{Sc}, the multimeric misfolded conformer PrP^{Sc} banded about half-way down the gradient. Discontinuous iodixanol gradients have also been used as part of the procedure in the purification of prion fibrils from mouse brain [38]. Sarkosyl-containing 10–54% (w/v) iodixanol gradients (200,000 g for 1 h) have also been used to determine the size of aggregates containing prion protein [39].

An extensive study of the quaternary structure of prion proteins was carried out by Laferrière et al [40] using both sedimentation velocity and equilibrium iodixanol density gradients. Mouse brain homogenates were solubilized in a buffered saline containing EDTA, DTT, dodecyl- β -D-maltoside and *N*-lauryl sarcosine.

- ◆ For sedimentation velocity analysis the samples were layered on top of a continuous 10–25% iodixanol gradient and centrifuged at 285 000 g for 45 min. In some studies the time was increased to 90 min; this provided additional resolution of the proteins. Rapidly- and slowly-sedimenting ovine strains were identified
- ◆ For equilibrium gradients the sample was adjusted to 40% iodixanol and made part of a discontinuous gradient spanning the range 10–60% iodixanol and centrifuged at 115,000 g for 17 hours.
- ◆ Similar sedimentation velocity gradients (5–25% iodixanol) were used by Coleman et al [41] in the identification of small soluble protease-sensitive prion proteins as the major source of infectivity in mouse brain.
- ◆ Wenborn et al [42] were able to purify prion protein from brain homogenates that was devoid of ferritin and other contaminating proteins by a combination of filtration and a series of centrifugations through 17.5% (w/v) iodixanol at 16,100 g.
- ◆ Herrmann et al [43] used 10–40% and 2–20% iodixanol gradients (52,000 rpm for 90 min) to show an impoverishment of the larger molecular weight prion proteins (from mouse brain) during drug treatment.
- ◆ Sarkosyl-solubilized brain material has been analyzed [44] in a 7–28% OptiPrep™ gradient (18,000 g for 20 min) and solubility assays were carried out in 15% OptiPrep™ at 18,000 g for 30 min [45, 46].
- ◆ Most recently, sedimentation velocity analysis has been carried out on 150 μ l of a sarkosyl-solubilized fractions loaded on to 4.8 ml 10–30% iodixanol gradients centrifuged at 285,000 g for 45 min [47].

9. Association of *Drosophila* proteins with lipophorin

This methodology was first reported by Eugster et al [48]: a homogenate of larvae was centrifuged at 120,000 g and the supernatant adjusted to adjusted to 50% (w/v) iodixanol and layered beneath gradients of 10%–20% to 35%–45% (w/v) iodixanol. They were centrifuged at 4°C, at 285,000 g, for 16 h in a swinging-bucket rotor (2.2 ml tube volume). There is no obvious reason why the separation cannot be scaled up. Proteins not bound to the lipophorin lipoprotein remained at the bottom of the gradient, while lipophorin-bound proteins floated into the top half of the gradient. The method was used to study the interaction of the lipophorin with various morphogens via heparan sulphate. Similar gradients for studying lipophorin interactions have been used by Palm et al [49, 50].

10. Hepatitis and Herpes virus proteins

To isolate hepatitis B surface antigen (sHBsAg) particles from *Leishmania* cells the latter were sonicated in a detergent-containing buffer and clarified by centrifugation at 8000 rpm, for 35 min. After particle formation (16–24 h) the supernatant was layered on top of a gradient of 6%, 12%, 18%, 24% and 30% (v/v) OptiPrep and centrifuged at 27,000 rpm for 16 h at 4 °C [51]. Hepatitis C subunits E1 and E2 have been analyzed on either sedimentation velocity or buoyant density gradients of 10–50%

iodixanol, centrifuged at 160,000 g_{av} for 16 h [52]. Herpes virus 8 interleukin 6 protein complexes were analyzed in gradients of 5, 10, 20, 30, 40 and 50% (w/v) iodixanol by centrifugation at 77,000 g_{av} for 3 h [53].

11. Other applications

Iodixanol gradients have also been used to study the size of complexes formed by the Circadian clock PERIOD complex [54] and using a small volume step gradient of 20%, 15%, 10%, 5%, centrifuged for 18 h at approx 150,000 g_{av} 4°C, Lin et al [55] were able to demonstrate an association of VASP and GST-WASP proteins from *Dictyostelium*. The gradient would be continuous within about 3 h of the centrifugation.

Some more recent publications include the following studies:

Alzheimers disease protein [56]
 β -Amyloid peptide – monomer separation [57]
Francisella tularensis secretion system structure [58]
Nicotiana benthamiana, protein expression in [59]
 ORF2 capsid protein [60]
 Prion proteins (intermolecular cross-linking) [61]
 (assembly in Parkinson’s disease) [62]
 (hydrophobic regions) [63]

11. References

1. Marcum, J.M. and Borisy, G.G. (1978) *Sedimentation velocity analyses of the effect of hydrostatic pressure on the 30 S microtubule protein oligomer* J. Biol. Chem., 253, 2852-2857 J. Biol. Chem., 253, 2852-2857
2. Hauge, J.G. (1971) *Pressure-induced dissociation of ribosomes during ultracentrifugation* FEBS Lett., 17, 168-172
3. Timasheff, S.N. (1993) *The control of protein stability and association by weak interactions with water: how do solvents affect these processes* Annu. Rev. Biophys. Biomol. Struct., 22 67-97
4. Rashid, D.J., Bononi, J., Tripet, B.P., Hodges, R.S. and Pierce, D.W. (2005) *Monomeric and dimeric states exhibited by the kinesin-related motor protein KIF1A* J. Pept. Res., 65, 538-549
5. Moir, R.D., Romano, D.M., Laurans, M.H., Huang, X., Bush, A.I., Smith, J.D. and Tanzi, R.E. (1999) *Differential effects of apolipoprotein E isoforms on metal-induced aggregation of A β using physiological concentrations* Biochemistry, 38, 4595-4603
6. Ward, R.V., Jennings, K.H., Jepras, R., Neville, W., Owen, D.E., Hawkins, J., Christie, G., Davis, J.B., George, A., Karran, E.H. and Howlett, D.R. (2000) *Fractionation and characterization of oligomeric, protofibrillar and fibrillar forms of β -amyloid peptide* Biochem. J., 348, 137-144
7. Rzepecki, P., Nagel-Steger, L., Feuerstein, S., Linne, U., Molt, O., Zadnart, R., Aschermann, K., Wehner, M., Schrader, T. and Riesner, D. (2005) *Prevention of Alzheimer’s disease-associated A β aggregation by rationally designed non-peptidic β -sheet ligands* J. Biol. Chem., 279, 47497-47505
8. Lockhart, A., Ye, L., Judd, D.B., Merritt, A.T., Lowe, P.N., Morgenstern, J.L., Hong, G., Gee, A.D. and Brown, J. (2005) *Evidence for the presence of three distinct binding sites for the thioflavin T class of Alzheimer’s disease PET imaging agents on β -amyloid peptide fibrils* J. Biol. Chem., 280, 7677-7684
9. Funke, S.A., van Groen, T., Kadish, I., Bartnik, D., Nagel-Steger, L., Brener, O., Sehl, T., Batra-Safferling, R., Moriscot, C., Schoehn, G., et al. (2010) *Oral treatment with the D-enantiomeric peptide D3 improves the pathology and behavior of Alzheimer’s disease transgenic mice* ACS Chem. Neurosci., 1, 639–648
10. Funke, S.A., Liu, H., Sehl, T., Bartnik, D., Brener, O., Nagel-Steger, L., Wiesehan, K. and Willbold, D. (2012) *Identification and characterization of an A β oligomer precipitating peptide that may be useful to explore gene therapeutic approaches to Alzheimer disease* Rejuvenation Res., 15, 144-147
11. Brener, O., Dunkelmann, T., Gremer, L., van Groen, T., Mirecka, E.A., Kadish, I., Willuweit, A., Kutzsche, J., Jürgens, D. et al (2015) *QIAD assay for quantitating a compound’s efficacy in elimination of toxic A β oligomers* Sci. Rep., 5: 13222
12. Rudolph, S., Klein, A.N., Tusche, M., Schlosser, C., Elfgen, A., Brener, O., Teunissen, C., Gremer, L., Funke, S.A., Kutzsche, J. and Willbold, D. (2016) *Competitive mirror image phage display derived peptide modulates amyloid beta aggregation and toxicity* PLoS One, 11: e0147470
13. Klein, A.N., Ziehm, T., Tusche, M., Buitenhuis, J., Bartnik, D., Boeddrich, A., Wiglenda T., Wanker, E., Funke, S.A. et al (2016) *Optimization of the all-D peptide D3 for A β oligomer elimination* PloS One, 11, e015035

14. Klein, A.N., Ziehm, T., van Groen, T., Kadish, I., Elfgen, A., Tusche, M., Thomaier, M., Reiss, K. et al (2017) *Optimization of D-peptides for A β monomer binding specificity enhances their potential to eliminate toxic A β oligomers* ACS Chem. Neurosci., 8, 1889–1900
15. Ziehm, T., Brener, O., van Groen, T., Kadish, I., Frenzel, D., Tusche, M., Kutzsche, J., Reiß, K., Gremer, L., Nagel-Steger, L. and Willbold, D. (2016) *Increase of positive net charge and conformational rigidity enhances the efficacy of D-enantiomeric peptides designed to eliminate cytotoxic A β species* ACS Chem. Neurosci., 7, 1088-1096
16. Van Groen, T., Schemmert, S., Brener, O., Gremer, L., Ziehm, T., Tusche, M. Nagel-Steger, L., Kadish, I., Schartmann, E. et al (2017) *The A β oligomer eliminating D-enantiomeric peptide RD2 improves cognition without changing plaque pathology* Sci. Rep., 7: 16275
17. Thomaier, M., Gremer, L., Dammers, C., Fabig, J., Neudecker, P. and Willbold, D. (2016) *High-affinity binding of monomeric but not oligomeric amyloid- β to ganglioside GM1 containing nanodiscs* Biochemistry 2016, 55, 6662–6672
18. Sehlin, D., Englund, H., Simu, B., Karlsson, M., Ingelsson, M., Nikolajeff, F., Lannfelt, L. and Pettersson, F.E. (2012) *Large aggregates are the major soluble Ab species in AD brain fractionated with density gradient ultracentrifugation* PLoS One, 7: e32014
19. Stöhr, J., Watts, J.C., Mensinger, Z.L., Oehler, A., Grillo, S.K., DeArmond, S.J., Prusiner, S.B. and Giles, K. (2012) *Purified and synthetic Alzheimer's amyloid beta (A β) prions* Proc. Natl. Acad. Sci. USA, 109, 11025-11030
20. Yeh, F.L., Wang, Y., Tom, I., Gonzalez, L.C. and Sheng, M. (2016) *TREM2 binds to apolipoproteins, including APOE and CLU/APOJ, and thereby facilitates uptake of amyloid-Beta by microglia* Neuron 91, 328–340
21. Streich, C., Akkari, L., Decker, C., Bormann, J., Rehbock, C., Muller-Schiffmann, A., Niemeyer, F.C., Nagel-Steger, L. et al (2016) *Characterizing the effect of multivalent conjugates composed of A β -specific ligands and metal nanoparticles on neurotoxic fibrillar aggregation* ACS Nano, 10, 7582-7597
22. Levy, M., Porat, Y., Macharach, E., Shalev, D.E. and Gazit, E. (2008) *Phenolsulfophthalein but not phenolphthalein inhibits amyloid fibril formation: implications for the modulation of amyloid self-assembly* Biochemistry, 47, 5896-5904
23. Frenzel, D., Glück, J.M., Brener, O., Oesterhelt, F., Nagel-Steger, L. and Willbold, D. (2014) *Immobilization of homogeneous monomeric, oligomeric and fibrillar A β species for reliable SPR measurements* PLoS One, 9: e89490
24. Khlistunova, I., Biernat, J., Wang, Y., Pickhardt, M., von Bergen, M., Gazova, Z., Mandelkow, E. and Mandelkow, E-M. (2006) *Inducible expression of Tau repeat domain in cell models of tauopathy* J. Biol. Chem., 281, 1205-1214
25. Mocanu, M-M., Nissen, A., Eckermann, K., Khlistunova, I., Biernat, J., Drexler, D., Petrova, O., Schönin, K., Bujard, H., Mandelkow, E., Zhou, L., Rune, G. and Mandelkow, E-M. (2008) *The potential for β -structure in the repeat domain of Tau protein determines aggregation, synaptic decay, neuronal loss and coassembly with endogenous Tau in inducible mouse models of Tauopathy* J. Neurosci., 16, 737-748
26. Li, T. and Paudel, H.K. (2016) *14-3-3 ζ Mediates tau aggregation in human neuroblastoma M17 cells* PLoS One, 11: e0160635
27. Shutova, M., Yang, C., Vasiliev, J.M. and Svitkina, T. (2012) *Functions of nonmuscle myosin II in assembly of the cellular contractile system* PLoS One, 7: e40814
28. Lehtimäki, J.I., Fenix, A.M., Kotila, T.M., Balistreri, G., Paavolainen, L., Varjosalo, M., Burnette, D.T. and Lappalainen, P. (2017) *UNC-45a promotes myosin folding and stress fiber assembly* J. Cell Biol., 216, 4053–4072
29. Lee, H-J. and Lee, S-J. (2002) *Characterization of cytoplasmic α -synuclein aggregates* J. Biol. Chem., 277, 48976-48983
30. MacCormick, M., Modersheim, T., van der Salm, L.W M., Moore, A., Pryor, S.C., McCaffrey, G. and Grimes, M.L. (2005) *Distinct signaling particles containing ERK/MEK and B-Raf in PC12 cells* Biochem. J., 387, 155-164
31. Betschinger, J., Eisenhaber, F. and Knoblich, J.A. (2005) *Phosphorylation-induced autoinhibition regulates the cytoskeletal protein lethal (2) giant larvae* Curr. Biol., 15, 276-282
32. Majoul, I.V., Bastiaens, P.I.H. and Soling H-D (1996) *Transport of an external Lys-Asp-Glu-Leu (KDEL) protein from the plasma membrane to the endoplasmic reticulum: studies with cholera toxin in Vero cells* J. Cell Biol., 133, 777-789
33. Chen, Q., Peto, C.A., Shelton, G.D., Mizisin, A., Sawchenko, P.E. and Schubert, D. (2009) *Loss of modifier of cell adhesion reveals a pathway leading to axonal degeneration* J. Neurosci., 29, 118-130

34. Minamide, L.S., Maiti, S., Boyle, J.A., Davis, R.C., Coppinger, J.A., Bao, Y., Huang, T.Y., Yates, J., Bokoch, G.M. and Bamburg, J.R. (2010) *Isolation and characterization of cytoplasmic cofilin-actin rods* J. Biol. Chem., 285, 5450-5460
35. Rademacher, S., Verheijen, B.M., Hensel, N., Peters, M., Bora, G., Brandes, G., de Sa, V.R., Heidrich, N., Fischer, S. et al (2017) *Metalloprotease-mediated cleavage of PlexinD1 and its sequestration to actin rods in the motoneuron disease spinal muscular atrophy (SMA)* Hum. Mol. Genet., 26, 3946–3959
36. Lin, W-H., Nelson, S.E., Hollingsworth, R.J. and Chung, C.Y. (2010) *Functional roles of VASP phosphorylation in the regulation of chemotaxis and osmotic stress response* Cytoskeleton, 67, 259–27
37. Tixador, P., Herzog, L., Reine, F., Jaumain, E., Chapuis, J., Le Dur, A., Laude, H. and Béringue, V. (2010) *The physical relationship between infectivity and prion protein aggregates is strain-dependent* PLoS Pathogens, 6: e1000859
38. Sim, V.L. and Caughey, B. (2009) *Ultrastructures and strain comparison of under-glycosylated scrapie prion fibrils* Neurobiol. Aging, 30, 2031–2042
39. Sigurdson, C.J., Joshi-Barr, S., Bett, C., Winson, O., Manco, G., Schwarz, P., Rüllicke, T., Nilsson, P.R., Margalith, I., et al (2011) *Spongiform encephalopathy in transgenic mice expressing a point mutation in the $\beta 2$ - $\alpha 2$ loop of the prion protein* J. Neurosci., 31, 13840–13847
40. Laferrière, F., Tixador, P., Moudjou, M., Chapuis, J., Sibille, P., Herzog, L., Reine, F., Jaumain, E., Laude, H., Rezaei, H. and Béringue, V. (2013) *Quaternary structure of pathological prion protein as a determining factor of strain-specific prion replication dynamics* PLoS Pathog., 9: e1003702
41. Coleman, B.M., Harrison, C.F., Guo, B., Masters, C.L., Barnham, K.J., Lawson, V.A. and Hill, A.F. (2014) *Pathogenic mutations within the hydrophobic domain of the prion protein lead to the formation of protease-sensitive prion species with increased lethality* J. Virol., 88, 2690-2703
42. Wenborn, A., Terry, C., Gros, N., Joiner, S., D'Castro, L., Panico, S., Sells, J., Cronier, S., Linehan, J.M., Brandner, S., Saibil, H.R., Collinge, J. and Wadsworth, J.D.F. (2015) *A novel and rapid method for obtaining high titre intact prion strains from mammalian brain* Sci. Rep. 5: 10062
43. Herrmann, U.S., Schütz, A.K., Shirani, H., Huang, D., Saban, D., Nuvolone, M., Li, B., Ballmer, B., Åslund, A.K.O. et al (2015) *Structure-based drug design identifies polythiophenes as antiprion compounds* Sci. Transl. Med., 7, 299ra123
44. Leske, H., Hornemann, S., Herrmann, U.S., Zhu, C., Dametto, P., Li, B., Laferriere, F., Polymenidou, M., Pelczar, P. et al (2017) *Protease resistance of infectious prions is suppressed by removal of a single atom in the cellular prion protein* PLoS One, 12: e0170503
45. Aguilar-Calvo, P., Xiao, X., Bett, C., Eraña, H., Soldau, K., Castilla, J., Nilsson, K.P.R., Surewicz, W.K. and Sigurdson, C.J. (2017) *Post-translational modifications in PrP expand the conformational diversity of prions in vivo* Sci. Rep., 7: 43295
46. Bett, C., Lawrence, J., Kurt, T.D., Orru, C., Aguilar-Calvo, P., Kincaid, A.E., Surewicz, W.K., Caughey, B., Wu, C. and Sigurdson, C.J. (2017) *Enhanced neuroinvasion by smaller, soluble prions* Acta Neuropath. Comm., 5: 32
47. Igel-Egalon, A., Moudjou, M., Martin, D., Busley, A., KnaEpple, T., Herzog, L., Reine, F. et al (2017) *Reversible unfolding of infectious prion assemblies reveals the existence of an oligomeric elementary brick* PLoS Pathog., 13: e1006557
48. Eugster, C., Panáková, D., Mahmoud, A. and Eaton, S. (2007) *Lipoprotein-heparan sulfate interactions in the Hh pathway* Devel. Cell, 13, 57-71
49. Palm, W., Sampaio, J.L., Brankatschk, M., Carvalho, M., Mahmoud, A., Shevchenko, A. and Eaton, S. (2012) *Lipoproteins in Drosophila melanogaster—assembly, function, and influence on tissue lipid composition* PLoS Genet., 8: e1002828
50. Palm, W., Swierczynska, M.M., Kumari, V., Ehrhart-Bornstein, M., Bornstein, S.R. and Eaton, S. (2013) *Secretion and signaling activities of lipoprotein-associated hedgehog and non-sterol-modified hedgehog in flies and mammals* PLoS Biol., 11: e1001505
51. Czarnota, A., Tyborowska, J., Peszyńska-Sularz, G., Gromadzka, B., Bieńkowska-Szewczyk, K. and Grzyb, K. (2016) *Immunogenicity of Leishmania-derived hepatitis B small surface antigen particles exposing highly conserved E2 epitope of hepatitis C virus* Microb. Cell Fact, 15: 62
52. Haddad, J.G., Rouillé, Y., Hanouille, X., Descamps, V., Hamze, M., Dabboussi, F., Baumert, T.F., Duverlie, G., Lavie, M. and Dubuisson, J. (2017) *Identification of novel functions for hepatitis C virus envelope glycoprotein E1 in virus entry and assembly* J. Virol., 91: e00048-17
53. Chen, D., Xiang, Q. and Nicholas, J. (2017) *Human herpesvirus 8 interleukin-6 interacts with calnexin cycle components and promotes protein folding* J. Virol., 91: e00965-17
54. Padmanabhan, K., Robles, M.S., Westerling, T. and Weitz, C.J. (2012) *Feedback regulation of transcriptional termination by the mammalian circadian clock PERIOD complex* Science, 337, 599-602

55. Lin, W-H., Nelson, S.E., Hollingsworth, R.J. and Chung, C.Y. (2010) *Functional roles of VASP phosphorylation in the regulation of chemotaxis and osmotic stress response* Cytoskeleton, 67, 259–271
56. Schemmert, S., Schartmann, E., Zafiu, K., Kass, B., Hartwig, S., Lehr, S., Bannach, O., Langen, K-J. et al (2019) *Aβ oligomer elimination restores cognition in transgenic Alzheimer's mice with full-blown pathology* Mol. Neurobiol., 56, 2211–2223
57. Mañucat-Tan, N.B., Saadipour, K., Wang, Y-J., Bobrovskaya, L. and Zhou, X-F. (2019) *Cellular trafficking of amyloid precursor protein in amyloidogenesis: physiological and pathological significance* Mol. Neurobiol., 56, 812–830
58. Ziveri, J., Chhuon, C., Jamet, A., Rytter, H., Prigent, G., Tros, F., Barel, M., Coureuil, M. et al (2019) *Critical role of a sheath phosphorylation site on the assembly and function of an atypical type VI secretion system* Mol. Cell. Proteom., 18, 2418–2432
59. Mbewana, S., Myers, A.E., Rybicki, E.P. (2019) *Chimaeric rift valley fever virus-like particle vaccine candidate production in Nicotiana benthamiana* Biotechnol. J., 14: 1800238
60. Ankavay, M., Montpellier, C., Sayed, I.M., Saliou, J-M., Wychowski, C., Saas, L., Duvet, S., Aliouat-Denis, C-M. et al (2019) *New insights into the ORF2 capsid protein, a key player of the hepatitis E virus lifecycle* Sci. Rep., 9: 6243
61. Teruya, K., Nishizawa, K., Oguma, A., Sakasegawa, Y., Kitamoto, T. and Doh-ura, K. (2019) *Intermolecular crosslinking of abnormal prion protein is efficiently induced by a primuline-sensitized photoreaction* Biochim. Biophys. Acta – Gen. Subj., 1863, 384–394
62. Igel-Egalon, A., Laferrière, F., Moudjou, M., Jan Bohl, J., Mezache, M., Knäpple, T., Herzog, L., Reine, F. et al (2019) *Early stage prion assembly involves two subpopulations with different quaternary structures and a secondary templating pathway* Commun. Biol., 2: 363
63. Abskharon, R., Wang, F., Wohlkonig, A., Ruan, J., Soror, S., Giachin, G., Pardon, E., Zou, W. et al (2019) *Structural evidence for the critical role of the prion protein hydrophobic region in forming an infectious prion* PLoS Pathog., 15: e1008139

OptiPrep™ Application Sheet M09; 8th edition, January 2020